

Now you can relax whenever you go out!

# MAKE YOUR DAY MORE REFRESHING!



The company pursues clean air in the nose, AirLab

**AirLab**

# 01 NOSECLEAN INTRODUCE

## Brand



**Brand name**  
NoseClean



**Brand Image**  
Clean and refreshing nose of human is embodied.

## Application of NoseClean use




The product is designed to help the user breathe clean air easily by filtering various pollutants from the air before it enters the nose, **including particulates, yellow dust, pollen, environmental dust, and car exhaust fumes**, via the fine hairs of a natural filter (hanji, or traditional Korean paper) and nano filter 864F.

## Example of wearing



## Order of NoseClean assembly



## How NoseClean differs from other products

- 'Inside nose filter mask' with 3-stage assembly configuration
- 'Environmentally-friendly filter' made with natural hanji and cotton
- 'Sanitary product' with disposable and replaceable single-use hanji 'gauze'
- Thin and transparent silicon creates a 'comfortable sense of wear' and 'minimal exposure'
- Addition of phytoncide to reinforce the **anti-bacterial effect** inside the nose
- **Refillable product** with reusable silicon tube and replaceable filter
- Individually customized products available in **various sizes** (product specifically designed for children, small, medium and large sizes)
- As the silicon tube is firmly attached to the nasal cavity, **only the filter is used when breathing.**
- In addition to the **primary filtration** using the fine hairs on the surface of the hanji and nano filter (filament: 1004EA), **secondary filtration** can also be provided using the moisture created inside the silicon (tube) due to the temperature difference in and outside the nose.

# 01 NOSECLEAN INTRODUCE

## Product Materials



**Main body**

LDPE, which is harmless to humans.



**Filter**

Cellulose (fabric) (hanji 58%) = natural filter  
Polyester [40%, filament 864 pieces]  
Cotton (2%) = natural filter



**Dual tube**

Silicon, which is harmless to humans.

## Product Composition

### Regular type

#### Regular type (3 pieces)

Children (Ø8mm), small (Ø9mm), medium (Ø10mm), large (Ø11mm)

- ① 3-stage assembly configuration, no flavor, odorless
- ② Refillable product with a replaceable natural filter (by reusing the silicon).
- ③ Comfortable wearing and minimal exposure are possible by means of thin and transparent silicon.
- ④ Primary filtration using the fine hairs on the surface of natural hanji and nano filter (filament: 864 pieces)

According to the result of a ultrafine particle collection efficiency test, the product blocked 70% for PM 2.5 or less (KIER).  
The result of an automobile exhaust gas test showed that the product blocked 20 to 30%.

#### Exchange filter regular type (12 sets)

- ① Natural filter:  
children and small (Ø5x6mm)  
medium (Ø6x7mm)  
large (Ø7x8mm)
- ② Material  
Cellulose (hanji content: 58%)  
= natural filter  
Polyester (40%) = filament 1004 pieces, ₩  
cotton (2%) = natural filter

### High-Quality type

#### High-Quality type (3 pieces)

- ① The product differs from conventional products because Ag (silver) has been added to the silicon to improve the anti-bacterial effect.
- ② Natural flavor and anti-bacterial action due to the addition of phytoncides to the natural filter.

#### High-quality exchange filter (12 sets)

- ① This product differs from conventional exchange filters in that the anti-bacterial effect is stronger and the addition of phytoncide to the filter creates a natural flavor.



## Product specifications



Model name	Small (Ø x H, mm)	Medium (Ø x H, mm)	Large (Ø x H, mm)	Distinction
① Regular type (NC-ST)	9 x 9	10 x 10	11 x 11	Silicon, which is harmless to humans
② High-quality type (NC-AG)	9 x 9	10 x 10	11 x 11	Silicon contains Ag (silver)
③ For children(regular) (NC-ST-K)	8 x 9			Silicon, which is harmless to humans
④ For children(high-quality type) (NC-AG-K)	8 x 9			Silicon contains Ag (silver)
⑤ Exchange filter(regular) (FT-ST)	5 x 6	6 x 7	7 x 8	Hanji (58%), Polyester (40%), Cotton(2%)/Odorless
⑥ Exchange filter(high-quality) (FT-AG)	5 x 6	6 x 7	7 x 8	Hanji (58%), Polyester (40%), Cotton(2%)/Addition of phytoncide



① 1 regular pack + 2 replaceable filter sets  
(Filter: Natural filter(odorless))  
small size, medium size



② Regular type 3 pieces  
(Filter: Natural filter(odorless))  
small, medium and large size



③ High-quality type 3 pieces  
(Filter: phytoncide, Silicon: Ag (silver))  
small, medium and large size



④ For children  
(Filter: Natural filter(odorless)),  
Silicon: transparent  
regular type



⑤ For children  
(Filter: Natural filter(odorless))  
Silicon: Ag (silver)  
high-quality type



⑥ Exchange filter (regular) 12 sets  
(Filter: Natural filter(odorless))  
kids&small, medium and large size




⑦ Exchange filter  
(Filter: phytoncide)  
small, medium and large size



# 02 USING NOSECLEAN

## How to use NoseClean

- 01 Form the silicon ring into a U-shape and insert deep into the nostrils.
- 02 Air is filtered via the fine hairs on the surface of the hanji and nano filter 1004F.
- 03 If the product falls from your nostrils, use a larger size.
- 04 If you feel that breathing quantity has been reduced, replace the filter.
- 05 Use from child use to adult use: small size is recommended for women, medium size for men.
- 06 When a child uses this product, the connecting ring may be visible from the outside.




## Precaution


- Before wearing this product, remove any foreign matter from inside your nose using a cotton swab.
- When eating food, remove the product and keep in its storage case.
- Wash the storage case once a day with water and completely dry before using.
- In the case of nasal discharge, first wipe your nose and then replace the filter before using.
- Incorrect wear results in poor air ventilation and other inconveniences. (Refer to the order of wearing)
- Before replacing the filter, wash the silicon with running water and remove only the moisture on the silicon before assembling.
- People who suffer from nasal stuffiness or a respiratory disorder must not use the product.
- In the case of an excessively runny nose, do not use the product.
- People with severe rhinitis must not use the product.
- NoseClean is not a medicine and, therefore, must not be taken orally.
- Do not go to sleep or enter water while wearing this product.
- Refrain from doing hard work or vigorous exercise.
- Do not use the product for any purpose other than the filtration of contaminated air.
- Keep the product out of reach of children.
- Do not use the product if the product for children does not fit into the nostrils of your child.

## How to use

### Product wearing procedures




### How to change the filter




## when purchasing the product

Check the size of your nostrils (actual interior measurement)

### Example of measurement method





※ Measure the size by placing a graduated ruler (K, S, M, L) under your nose.



	Size	Silicon diameter (Ø) x Height(H)
	For children (Kids)	8x9
	Small (S)	9x9
	Medium (M)	10x10
	Large (L)	11x11

unit: mm






① 1004F  
Top of 1004F Filter  
(1st filtration)

②  
Generation of  
humidity caused  
by temperature  
difference in the nose  
(2nd filtration)

#### Detailed drawing of Nose Clean 2 Filter



#### Test Report (for each filter)

Test Items	Test Result
No. of Filaments (EA):: KS K 0215:2012	216 strands



# 03 PRODUCT TEST REPORT

## PRODUCT TEST REPORT

**시험성적서**

1. 분 적 서 번 호 : KCT-02947  
2. 의 례 자 : ○ 공 격 일 : 2017년 11월 17일  
3. 시 행 기 간 : 2017년 11월 17일 ~ 2017년 11월 17일  
4. 시 행 장 소 : 서울특별시 강남구 테헤란로 157 (삼성전자빌딩 157층)  
5. 시 행 용 : 노즈클린 필터  
6. 시 행 품 명 : (1) KCT-02 1105-2-2015  
7. 시 행 결 과 : (1) 노즈클린 필터


시행품명	시험방법	시험결과	시험일자	시험장소	시험인원
노즈클린 필터 (KCT-02 1105-2-2015)	노즈클린 필터 시험	46.7	2017.11.17	서울특별시 강남구 테헤란로 157 (삼성전자빌딩 157층)	김민준, 김민준

2017년 11월 17일

한국건설생활환경시험연구원

### Nose Clean Filter Filter test (gas filtration)

〈KCL Test Result 2017.11.17〉



—■— Nose Clean Filter —●— KF 80 Filter

**시험 결과**

시험자 번호: KIER-16-0335호  
제1차시험 / (중7)

**항목 4. 필터효율 산출 방법**

○ 시험용 마스크를 착용한 피검자의 호흡량을 측정하고, 마스크 후단에서의 미세먼지 농도를 측정하고, 4회 측정된 후 평균값을 구하여 아래의 수식 (1)을 이용하여 피검자의 총량 필터효율 (%)을 산출함.

$$\eta_f = 1 - \frac{C_{in} - C_{out}}{C_{in}} \times 100, \% \quad (1)$$

○ 시험용 피검자용 마스크가 고장된 필터효율의 확인에 입자 손실이 확인되지 않거나, 호흡을 측정하기 위해 시험용 마스크의 비누제를 제거한 후, 필터효율 시험을 동일하게 4회 반복하여 수식 (1)을 이용하여 입자손실 효율 (%)을 계산함.

○ 필터효율 확인에 입자 손실이 보정된 필터효율 (%)을 아래의 수식 (2)을 이용하여 계산함.

$$\eta_b = \eta_f - \eta_{in} \quad (2)$$

○ 유량별 피검자용 마스크의 미세먼지 (PM2.5) 필터효율

No	유량 (L/min)	$\eta_f$	$\eta_{in}$	$\eta_b$
1	4.0	69.1	21.8	47.3
2	8.0	58.3	26.7	31.6

### Test Request : Filter dust collection efficiency

KIER (Korea Institute of Energy Research)

〈KIER Test Result : 2016. 06. 09〉

## 1. Test condition

### ① Test particles:

Particles generated by an Atomizer (9302, TSI, USA)

### ② Fine dust concentration :

Applicable only to aeromechanical diameter of 2.5um or fine dust concentration of PM2.5.

## 2. Test results

(fine dust PM 2.5 dust collection efficiency testing)

Flux 8.0, when breathing

**58.3% filtration**

Flux 4.0, when breathin

**69.1% filtration**



Front  
(nasal cavity)



Back  
(entrance)