




Tricon
New Series




Dark Gray




Light Gray




Hazel Brown




Maroon Brown



Turquoise



Gold





Intense

DLA 14.0mm / G.DIA 13.5mm / BC 8.6mm

Ice



DIA 14.5 mm | G.DIA 13.5 mm | BC 8.6 mm



Gray



Brown




Blue




Green

3078


DIA 14.5mm / G.DIA 13.9mm / BC 8.6mm




3078 Gray



3078 Brown



3078 Blue




3078 Green




Iricon Ravishing

— 1 day / Monthly —




Ash Gray

DIA : 14.5mm
G.DIA : 15.5mm
BC : 8.6mm




Wood Brown

DIA : 14.5mm
G.DIA : 15.5mm
BC : 8.6mm



Deep Blue

DIA : 14.5mm
G.DIA : 15.5mm
BC : 8.6mm



Jade Green

DIA : 14.5mm
G.DIA : 15.5mm
BC : 8.6mm



Iricon Lovesome Monthly




Monthly / DIA 14.2mm / G.DIA 13.0mm / BC 8.6mm / 2PCS/PACK



Lilian Brown

*Lovely eyes with a pink brown color
that permeates smoothly*

DIA 14.2mm | G.DIA 13.0mm | B.C 8.6mm



Light pupil



Medium pupil




Dark pupil



Pure Gray

*The perfect combination of
the soft blue and gray colors
that give you a refreshing look creates deep eyes*

DIA 14.2mm | G.DIA 13.0mm | B.C 8.6mm



Light pupil



Medium pupil



Dark pupil



Lens Specification

Material	HEMA
Packing	2PCS/PACK
Diameter	14.2mm(DIA) / 13.0mm(G.DIA)
Base Curve	8.6mm
Water Contents	42%
Power Range	$\pm 0.00 \sim -10.00$
Use Cycle	1month
Color	Lilian Brown / Pure Gray



EYEMADE



Vision Vibe



CLEAR2SEE

IRICON 1-DAY SERIES



DIA 14.5mm G.DIA 13.8mm

|

DIA 14.5mm G.DIA 13.5mm

EYEMADE

Lens Specification

Material	HEMA
Packing	10pcs/pack, 30pcs/pack
Diameter	14.5mm(DIA) / 13.8mm(G.DIA)
Base Curve	8.6mm
Water Contents	42%
Power Range	± 0.00 ~ -10.00
Use Cycle	1DAY
Color	London Gray,Tokyo Brown,Oslo Blue,Mexico Green, Paris Violet



Vision Vibe

Lens Specification

Material	HEMA
Packing	10pcs/pack, 30pcs/pack
Diameter	14.5mm(DIA) / 13.5mm(G.DIA)
Base Curve	8.6mm
Water Contents	42%
Power Range	± 0.00 ~ -10.00
Use Cycle	1DAY
Color	Luna Gray,Venus Brown,Mercury Blue,Saturn Green,Mars Violet



CLEAR2SEE


Lens Specification

Material	HEMA
Packing	10pcs/pack, 30pcs/pack
Diameter	14.0mm(DIA)
Base Curve	8.6mm
Water Contents	42%
Power Range	± 0.00 ~ -10.00
Use Cycle	1DAY
Color	Blue tint



1DAY IRICON EYEMADE

DIA 14.5mm G.DIA 13.8mm BC 8.6mm



POINT 1

Gentle color that makes your eyes look naturally


POINT 2

Two-tone color lens with a gradient design that makes you fall in love with

POINT 3

Clear ring circle line to create a lively eyelook

EYEMADE COLOR




London Gray



Tokyo Brown



Oslo Blue



Mexico Green




Paris Violet



Iricon Vision Vibe

DIA 14.5mm G.DIA 13.5mm BC 8.6mm







Luna Gray




Venus Brown



Mercury Blue



Saturn Green



Mars Violet











IRICON CLEAR2SEE

ONE DAY




IRICON COSPRÉ LENS




Clear eyes that look bigger naturally

Iricon lens allows you to choose a product according to the color and size of your iris, and its natural ring circle line makes your eyes look bigger and clearer. Experience different designs with different patterns and colors !



G.DIA


[13.4]



2015

Clear & lively


[13.5]



2006

Falling in love with


[13.6]



2039

Delicate but clear


[13.8]



2050

Natural & lively

[14.0]



2023

Calm & clear


[13.0]



3046

Colorful & unique


[13.2]



3089

Dreamy eye


[13.4]



3064

Gentle but brilliant


[13.7]



3005

Mysterious

[13.9]



3081

Adorable & bubbly

1 - tone color



1001A



1001B




1001C



1001D



1001E




1002A



1002B



1002C




1002D



1002E

2 - tone color




2001A




2001B




2001C




2001D



2001E




2002A




2002B




2002C




2002D




2002E




2003A




2003B




2003C



2003D



2003E



2004A



2004B




2004C



2004D



2004E



2005A



2005B




2005C



2005D




2005E




2006A



2006B




2006C



2006D



2006E



2007A



2007B



2007C



2007D



2007E



2008A



2008B



2008C



2008D



2008E




2009A



2009B



2009C




2009D



2009E



2010A



2010B




2010C




2010D




2010E




2011A




2011B



2011C




2011D



2011E




2012A



2012B




2012C



2012D



2012E



2013A



2013B




2013C




2013D




2013E




2014A




2014B



2014C



2014D




2014E




2015A




2015B



2015C




2015D




2015E




2016A



2016B



2016C



2016D



2016E



2017A



2017B




2017C



2017D




2017E




2018A




2018B




2018C



2018D



2018E




2019A




2019B




2019C




2019D




2019E




2020A




2020B



2020C




2020D



2020E


3 - tone color




3001A



3001B



3001C



3001D



3001E



3002A



3002B




3002C



3002D




3002E




3003A




3003B




3003C



3003D



3003E



3004A



3004B



3004C



3004D



3004E



3005A



3005B



3005C



3005D



3005E



3006A



3006B




3006C



3006D



3006E






3017A




3017B



3017C



3017D




3017E




3018A



3018B




3018C




3018D




3018E




3020A




3020B




3020C




3020D




3020E




3021A




3021B




3021C



3021D



3021E



3022A




3022B



3022C



3022D



3022E

4 - tone color




4001A



4001B



4001C



4001D



4001E

Fancy Lenses



9001



9002



9003




9004



9005


Icon Color Chart

Skin tone




COOL TONE


LIGHT




**EYEMADY
LONDON GRAY**




**LOVESOME MONTHLY
PURE GRAY**




**3056
BLUE**




**LOVESOME MONTHLY
LILIAN BROWN**



**3010
BLUE**



**VISION VIBE
LUNA GRAY**



**2006
GRAY**

SUMMER

WINTER

DEEP

Iricon Color Chart

Skin tone



WARM TONE

LIGHT



**3039
BROWN**



**3089
BROWN**




**EYEMADY
TOKYO BROWN**



**VISION VIBE
VENUS BROWN**

SPRING

AUTUMN



**3016
GRAY**



**3012
BROWN**

DEEP

Iricon

DIA : 14.5mm



Iricon

DIA : 14.2mm




Iricon

DIA : 14.0mm



Iricon color



1 DAY



- Suitable for light & dark eyes
- All color & size available

1 MONTH



- Suitable for light & dark eyes
- All color & size available

3 MONTHS




- Suitable for light & dark eyes
- All color & size available

3-Layer Process

Safety is important so we insist colorants do not come into contact with the cornea or eyelid.

Pigments embedded in the middle layer.



Optimum Water Content

42% water content keeps long-lasting moisture and high breathability. The lens is soft and locks in moisture and more comfortable.



Special Structure

Space is reserved between eyeball and lens edge, ensuring that tears can go through smoothly bringing in more oxygen. Our own special design ensures accurate center orientation so that wearers get better vision.



UV Blocking

Iricon Color blocks UV radiation and protects the eyes from frontal exposure and oblique, scattered UV rays from the sides.



Iricon clear

Level up your eye health with Iricon-Clear Lens

1 DAY



- 8.6 / 14.0, 14.2 Available
- Light-Blue Tinted

WEEKLY



- 8.6 / 14.0, 14.2 Available
- Light-Blue Tinted

MONTHLY



- 8.6 / 14.0, 14.2 Available
- Light-Blue Tinted

High Oxygen Permeability

Contact lens is an artificial device whose front surface substitutes the anterior surface of the cornea.

So it is very important to keep your eye with fresh oxygen.

Freely oxygen permeable eliminate potential damage to your cornea and enrich oxygen permeability.

The Best HEMA Material

Iricon-Clear can provide comfort because of material softness and keep maintaining hydration and wettability.

It gives you better vision because of its softness.

Special Design

To prevent and minimize deposits, Iricon-Clear use our own special design.

It can help tear and water circulation in eyes, and prevent or reduce the risk of ocular infection.

UV Protection

We provide you with protection against damaging rays. Iricon-Clear can fully protect the eyes from harmful light, and it is enough to keep your eyes safe from sun damage without other precaution.



IRICON ♥ LENS

SPECIAL CUSTOMIZATION

G.DIA?

**You can adjust
the graphic size of
the design you want to
fit your eyes**

Size?

**14.0, 14.2, 14.5 DIA
available in all design**

Color?

**Brown, gray, blue, green,
and violet are used to create the colors I want**

Production?

**New products in good condition
by producing on-demand**





Fascinate Design



Direct Manufacturing



Assured Quality



Best Material



Reasonable MOQ



Best Price



Fast Delivery



Service Support



Customization

KWIN Co., Ltd.

Factory & Headquarter

59, Sandanoryong-Gil, Samgi-Myeon, Iksan-Si,
Jeonbuk-State, Republic of Korea

Seoul R&D Center (Overseas Inquiry)

5th Floor, 11, Yeonhui-Ro, Mapo-Gu, Seoul

E-mail: kwin101@kwin.or.kr

www.kwinglobal.com





Iricon



Contact lens blossom, Iricon!
The world most fascinate design contact lens,
and well-made contact lens